



# CANADIAN ANIMAL BLOOD BANK

2055 Notre Dame Ave. Room AB 71

Winnipeg, MB R3H 0J9

(204) 632-2586

## CROSSMATCH TECHNIQUE

A crossmatch is performed to detect serological incompatibility.

### **DOGS**

Crossmatching is recommended before transfusion, although dogs usually lack naturally occurring alloantibodies. All dogs that have received red cells more than four days previously must be crossmatched with a fresh sample of serum before receiving any more red cells.

### **CATS**

Cats have naturally occurring alloantibodies and a crossmatch should **always** be performed prior to transfusion.

### **Equipment and Reagents**

5 glass tubes<sup>1</sup>

Test tube rack(s)

Pasteur pipettes (plastic or glass)

Incubator or waterbath at 37° C<sup>3</sup>

Centrifuge

Phosphate Buffered Saline (PBS)<sup>2</sup>

Scissors

Thermometer

### **Procedure**

1. Prepare 37°C waterbath. If using a styrene box, this may take 30 minutes to achieve equilibrium to hold the temperature for the 30 minute incubation time.
2. Collect blood from the RECIPIENT in a tube containing no anticoagulant or serum separator.
3. Centrifuge the specimen at 1000 g for at least 5 minutes to separate the serum from the red cells.
4. Remove the serum with a pipette into a clean, labelled test tube.
5. Remove some of the free red cells at the bottom of the original blood collection tube into a clean labelled test tube.
6. Prepare DONOR red cells to be crossmatched as follows:
  - a) Cut off a segment at the heat seal from a unit of red cells concentrate and place into a test tube labelled with the unit number.
  - b) Cut off one end of the red cell segment. Invert the segment so that the open end is in the test tube.
  - c) Cut off the other end of the segment. The red cells should empty into the test tube.

---

<sup>1</sup> 10mm x 75mm glass tubes (or equivalent) – enough tubes for 1 crossmatch, add 2 more tubes for each additional crossmatch unit

<sup>2</sup> Although 0.9% saline can be used, PBS reacts better with canine and feline red cells.

<sup>3</sup> A waterbath can be made with a styrene container with a ceramic mug (heat sink), will hold a test tube rack and filled with 2 - 3 inches of tap water at 37° C.

- d) Wash the recipient and donor red cells with saline.
- Fill the test tube with saline.
  - Centrifuge tube at 1000 g for 1 minute.
  - Decant the supernatant.
  - Repeat 2 more times.
  - If the supernatant shows haemolysis, continue washing the cells until the supernatant is clear and colourless.
- e) Make a 5% suspension of the red cells in the same tube.  
(1 drop of red cells + 20 drops saline)

7. Prepare tubes as follows for EACH unit to be crossmatched:

Label tube with:	<u>CANINE</u> 37° C Phase
Auto Control	4 drops <b>RECIPIENT</b> serum + 1 drop of 5% <b>RECIPIENT</b> red cell suspension (see 5.e)
Unit Number to be crossmatched	4 drops <b>RECIPIENT</b> serum + 1 drop of 5% <b>DONOR</b> red cell suspension

Label tube with:	<u>FELINE ONLY</u> Room Temp Phase
Auto Control	4 drops <b>RECIPIENT</b> serum + 1 drop <b>RECIPIENT</b> red cells suspension
Donor to be crossmatched	4 drops <b>RECIPIENT</b> serum + 1 drop <b>DONOR</b> red cell suspension
Minor Crossmatch of Donor Unit	4 drops <b>DONOR</b> serum + 1 drop <b>RECIPIENT</b> red cell suspension
Donor Control	4 drops <b>DONOR</b> serum + 1 drop <b>DONOR</b> red cell suspension

8. Mix all tubes gently.
9. Incubate the CANINE set of tubes at 37° C for 30 minutes. (Check the temperature of the incubator).

<ul style="list-style-type: none"> <li>• <b>PROCEED DIRECTLY TO STEP 10 FOR <u>FELINE</u> TUBES!</b></li> </ul>
---

10. Centrifuge all tubes for 15 seconds at 1000 g.
11. Examine the supernatant for haemolysis.
12. **GENTLY** suspend the cell button by rocking the tube and examine for macroscopic agglutination.

### Interpretation

**Compatibility** is indicated when there is no haemolysis or agglutination in any test tube, either macroscopically or microscopically. Test result is NEGATIVE.

**Incompatibility** is indicated by any haemolysis and/or agglutination in the red cells crossmatch tubes but not in the control tubes. A new donor must be selected. Test result is POSITIVE. If the control cells are weakly positive and the test sample is a strong positive, the crossmatch is **incompatible**.

If the control cells and the test cells are equally positive, **no conclusions as to compatibility can be drawn**.

### Notes

A compatible crossmatch does not prevent sensitization or a delayed transfusion reaction. It simply indicates that at the present time there are no significant antibodies against the red cells.